

Primer pairs used for each *M. pneumoniae* MLST locus

Name		Primer sequence (5'-3')	Amplicon (bp)	Location of MLST locus in CDS	Tm (°C)
<i>ppa</i>	F	CGCTGACCAAGCCTTTCTAC	256	192-440	60
	R	CACTCCAAACTTTGCACTCCC			
<i>pgm</i>	F	AGCACCTTGCACGATGAAGA	1072	456-1652	60
	R	CCTGCGCCTTCGTTAATTGG			
<i>gyrB</i>	F	TTGTCCCGGACTTTACCGTG	429	524-952	60
	R	TGTTTTCGACAGCAAAGCGG			
<i>gmk</i>	F	GAGCGGTGTTGGCAAAAGTA	394	189-582	60
	R	TGCATCCTCGTCATTACGCTT			
<i>glyA</i>	F	CAGAGAACTATGTGAGTAGGGACA	676	74-749	60
	R	TGACAACCCGGAAAGACACC			
<i>atpA</i>	F	GTCGCTGATGGCATTGCTAAG	796	100-895	60
	R	CCAGTAAACGCGAGTGCAAG			
<i>arcC</i>	F	CCCCATCAAGCCGTGTA	570	304-873	60
	R	TTGGGCAATAATGGCCGTCT			
<i>adk</i>	F	GTAGCCAACACCACCGGATT	473	70-542	60
	R	ACGGTGTCTTCGTAAAGCGT			